CUSTOM MOLDED EARPLUG - DECIBULLZ CUSTOM MOLDED EARPLUG, RED

Save Your Hearing! Custom Molded Ear Plugs

The Decibullz custom molded earplugs mold easily and fit the exact shape of your ear. This creates a perfect fitting earplug that is not painful, never falls out, and provides superior noise isolation. If you don't get the perfect fit the first time Decibullz are the only custom earplugs that are re-moldable. Weight 0.15 lbs
Rating: Meet ANSI S12.6 (ANSI s3.9) with independently tested NRR of 31
Included: Two Decibullz Thermoplastic Custom Molds, Three Sets of Triple Flange Tips S, M, L, One Set of Max Protect MEDIUM Foam Tips, Earplug Carrying Pouch with Decibullz Logo
Decibullz Earplugs are perfect for any situation where you need to protect your hearing, including shooting sports, DIY home improvement projects, and manufacturing. The earplugs have the highest rating of any custom earplug on the market with a Noise Reduction Rating of 31.

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Safety Instruction Guide for CUSTOM MOLDED EARPLUG DECIBULLZ

Introduction

Thank you for choosing the Decibullz Custom Molded Earplugs. These earplugs are designed to provide superior noise isolation and comfort. This guide contains important safety information and instructions to ensure safe and effective use of the product.

General Safety Guidelines

- Ensure that the earplugs are used only by individuals who are 12 years of age or older.
- Always inspect the earplugs before each use for any signs of damage or wear.
- Keep the earplugs out of reach of small children and pets to prevent choking hazards.
- Do not use the earplugs in situations where hearing is critical (e.g., swimming, cycling, or operating machinery) unless you are aware of your surroundings.
- If you experience any discomfort or pain while using the earplugs, remove them immediately and consult a healthcare professional if necessary.

Specific Safety Precautions for Use

- Always ensure a proper fit before use. If the earplugs do not fit correctly, remold them according to the instructions provided.
- Do not share earplugs with others to prevent ear infections or transmission of bacteria.
- Avoid using earplugs if you have an ear infection or a history of ear problems without consulting a healthcare professional.
- Do not use the earplugs while sleeping unless they are specifically designed for that purpose.
- Limit the duration of use to prevent discomfort or irritation in the ear canal.

Instructions for Installation and Usage

1. Molding the Earplugs:

- Prepare the Decibullz Thermoplastic Custom Molds as per the instructions included in the package.
- Heat the molds in hot water until they become pliable.
- Insert the molds into your ears and hold them in place until they cool and harden. This will create a custom fit.

2. Inserting the Earplugs:

- Ensure your hands are clean before handling the earplugs.
- Gently pull back on the outer ear to open the ear canal.
- Insert the molded earplug into the ear canal until it fits snugly.
- Repeat the process for the other ear.

3. Removing the Earplugs:

- Grasp the earplug and gently twist to break the seal.
- Pull the earplug out slowly and carefully to avoid discomfort.

4. Cleaning and Maintenance:

- Clean the earplugs regularly with mild soap and water.
- Allow them to dry completely before storing them in the provided carrying pouch.

Disposal Instructions

- Dispose of the earplugs responsibly. If they are damaged or no longer provide a proper fit, discard them in the regular waste.
- Do not dispose of the earplugs in water bodies or environments where they may cause harm to wildlife.

Contact Information for Further Support

For any safety inquiries, product concerns, or further assistance, please reach out to the appropriate customer support channels provided by the manufacturer.

By following these guidelines, you can ensure the safe use of your Decibullz Custom Molded Earplugs while enjoying superior noise isolation and comfort. Thank you for prioritizing your hearing protection.

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